

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(styrene-co-methacrylic acid)

Product Number: **678240**

Formula: $(C_8H_8)_x(C_4H_6O_2)_y$



TEST

Specification

Appearance (Color)	White
Appearance (Form)	Powder or Granules
Infrared spectrum	Conforms to Structure
Glass Transition by DSC	95.0 - 99.0 °C
Gel Permeation Chromatography	Corresponds to Requirements
Ratio	Confirmed
Molar Ratio 10:1	

Specification: PRD.1.ZQ5.1000098411

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.