


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Griseofulvin - from *Penicillium griseofulvum*, 97.0-102.0%

Product Number: **G4753**
CAS Number: 126-07-8
MDL: MFCD00082343
Formula: C₁₇H₁₇ClO₆
Formula Weight: 352.77 g/mol
Storage Temperature: -20 °C



TEST

Specification

| | |
|---------------------------|-----------------------------|
| Appearance (Color) | White to Off-White |
| Appearance (Form) | Powder |
| Solubility (Color) | Very Faint Yellow to Yellow |
| Solubility (Turbidity) | Clear |
| 50 mg/ml, DMF | |
| Infrared spectrum | Conforms to Structure |
| Assay | 97.0 - 102.0 % |
| Recommended Retest Period | ----- |
| 6 years | |

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.