

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Carbon nanotube, single-walled - 98% (Semiconducting)

Product Number: 750522
CAS Number: 308068-56-6

C

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Solid
May appear as a thick film powder or sheets	
SWNT	Confirmed
Semi-Conducting Purity	
Specification: Report Result	
Metallic Purity	
Specification: Report Result	

Specification: PRD.0.ZQ5.10000030070

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.