

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Collagen solution from bovine skin – BioReagent, suitable for cell culture, sterile-filtered

Product Number: C4243
CAS Number: 9007-34-5
MDL: MFCD00130825

Storage Temperature: 2 - 8 °C

TEST	Specification
Concentration Collagen	2.9 - 3.2 mg/ml
pH	1.9 - 2.1
Gel Time < or = 40 minutes	Conforms
Sterility No growth	Conforms
Endotoxin Assay By LAL, < or = 1.0 EU/mL	Conforms
Expiration Dated Product	-----

Specification: PRD.2.ZQ5.10000018350

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.