

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

HydroMatrix™ Peptide Cell Culture Scaffold – mixture, powder

Product Number:

A6982

Formula Weight:

1,526.47 g/mol

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Powder

Molecular Weight

Conforms

(by Mass Spec) Exhibits correct molecular weight

Purity (HPLC)

≥ 92 %

Specification: PRD.0.ZQ5.10000044683

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.