


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium 2-methoxypyrimidine-5-trifluoroborate - 97%

Product Number: 728810



Formula: C₅H₅BF₃KN₂O
Formula Weight: 216.01 g/mol

TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals
Proton NMR spectrum	Conforms to Structure
Fluorine NMR Spectra	Conforms to Structure
Carbon	26.8 - 28.8 %
Nitrogen	12.5 - 13.5 %
ICP Major Analysis	Confirmed
Confirms B and K Components	
Purity (HPLC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000052875

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.