

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Phosphodiesterase I from Bothrops atrox - Type V, crude dried venom,  $\geq 0.01$  unit/mg solid

Product Number:

**P4631**

CAS Number:

9000-89-9

MDL:

MFCD00131883

Storage Temperature:

-20 °C

**TEST****Specification**

units/mg solid

 $\geq 0.01$ 

Enzymatic Activity Unit Definition: One unit will hydrolyze 1.0 micromole of bis(p-nitrophenyl)-phosphate per minute at pH 8.8 at 37 deg C.

Specification: PRD.0.ZQ5.10000077530

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.