

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
MEH-PPV - average M_n 40,000-70,000

Product Number: **541443**
CAS Number: 138184-36-8
MDL: MFCD03093943



Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Conforms to Requirements
Red or Orange-Red	
Appearance (Form)	Powder or Fibers
Infrared spectrum	Conforms to Structure
Solubility (Turbidity)	Clear to Slightly Hazy
c = 0.02 g/l; Toluene	
Solubility (Color)	Conforms to Requirements
Bright Orange	
Molecular Number by GPC	40000 - 70000

Specification: PRD.2.ZQ5.10000020360

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.