


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N,N-Diethyl-1,1,2,3,3,3-hexafluoropropylamine

Product Number: 564990
CAS Number: 309-88-6
MDL: MFCD00054683
Formula: C₇H₁₁F₆N
Formula Weight: 223.16 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 90.0 %

Note: GC analysis includes Et₂NCF=CF₂CF₃ (two isomers), both of which are active fluorinating agents.

Specification: PRD.1.ZQ5.10000028694

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.