

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Formaldehyde solution - for molecular biology, 36.5-38% in H2 O

**Product Number:** F8775 CAS Number: 50-00-0

MDL: MFCD00003274

Formula: CH2O

Formula Weight: 30.03 g/mol



## TEST Specification

Color Test < 10 APHA
Suitability Suitable

Suitable For Use In Formaldehyde-Agarose Gel Electrophoresis of RNA

Recommended Retest Period

2 Years

Specification: PRD.0.ZQ5.10000026007

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.