

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

2 Years

dCas9-3XFLAG™-Biotin Protein

Product Number: DCAS9PROT

Storage Temperature: -20 ℃

TEST	Specification	
Vial Content (Overpack)	≤ 50.00 %	
Purity	≥ 90 %	
SDS-PAGE		
Activity	Conforms	
gRNA dependent gel shift assay		
Non-Specific DNase Activity	None Detected	
Non-Specific Nickase Activity	None Detected	
RNase Activity	None Detected	
Biotin Incorporation	Conforms	
Endotoxin Level	< 1.00 EU/ug	
Non-Specific DNase Activity	None Detected	
DBUFFER		
RNase Activity	None Detected	
DBUFFER		
Endotoxin Level	≤ 1.00 EU/ml	
DBUFFER		
Non-Specific DNase Activity	None Detected	
RSOLUTION		
RNase Activity	None Detected	
RSOLUTION		
Endotoxin Level	< 1.00 EU/ml	
RSOLUTION		
Recommended Retest Period		

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.



3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: dCas9-3XFLAG™-Biotin Protein

Product Number: DCAS9PROT

Storage Temperature: -20 ℃

TEST Specification

Specification: PRE.1.ZQ5.10000016967

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.