


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N,N'-Hexamethylene bis(acetamide) - 98%

Product Number: 224235
CAS Number: 3073-59-4
MDL: MFCD00008684
Formula: C₁₀H₂₀N₂O₂
Formula Weight: 200.28 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Crystalline Powder or Flakes and/or Lumps	Conforms to Requirements
Proton NMR Spectrum	Conforms to Structure
Carbon	58.5 - 61.5 %
Nitrogen	13.6 - 14.3 %
Purity (GC)	≥ 97.5 %

Specification: PRD.2.ZQ5.10000031367

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.