


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ammonium hydrogen difluoride - 99.999% trace metals basis

Product Number: 455830
CAS Number: 1341-49-7
MDL: MFCD00012018
Formula: H5F2N
Formula Weight: 57.04 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Crystals
X-Ray Diffraction	Conforms to Structure
Gravimetric Analysis % F with Lead Acetate	64.9 - 68.3 %
Purity 99.999% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.0.ZQ5.10000068713

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.