


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,2-Diphenylethanol – 99%

Product Number: 196568
CAS Number: 1883-32-5
MDL: MFCD00004729
Formula: C₁₄H₁₄O
Formula Weight: 198.26 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
GC (area %)

White
Powder
Conforms to Structure
≥ 98.5 %

Specification: PRD.0.ZQ5.10000013062

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.