

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

**Product Specification**

Product Name: Monoclonal Anti-β-Tubulin III (neuronal) antibody produced in mouse – ~1.0 mg/mL, clone 2G10, purified immunoglobulin, buffered aqueous solution  
Product Number: T8578

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.0 - 1.5 mg/ml
WB-Cell Line/Tissue Extract	PC-12 cells
WB Working Concentration	0.25 - 0.50 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000016379

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.