

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum(0)-1,3-divinyl-1,1,3,3-tetramethyldisiloxane complex solution - 0.1 M in poly(dimethylsiloxane), vinyl terminated


Product Number:

479527

MDL:

MFCD00151662

Formula:

C₈H₁₈O₂PtSi₂

TEST

Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

ICP Major Analysis

Confirmed

Confirms Pt and Si Components

Assay

≥ 1.8 %

Platinum

Specification: PRD.0.ZQ5.10000019962

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.