


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Aminobenzoic acid - ReagentPlus®, 99%

Product Number: 100536
CAS Number: 150-13-0
MDL: MFCD00007894
Formula: C₇H₇NO₂
Formula Weight: 137.14 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	98.5 - 101.5 %
Purity (HPLC)	≥ 98.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
c = 5%, 95% Ethanol	
Solubility (Color)	Colorless to Yellow

Specification: PRD.1.ZQ5.10000017833

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.