

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Monoclonal Anti-m-Calpain (Domain III/IV) antibody produced in mouse - clone 107-82, ascites fluid

Product Number: **C268**

MDL: MFCD04041214

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

ascites fluid containing 0.05% sodium azide

Working Dilution	Conforms
------------------	----------

By Direct Immunoblotting; 1:500

Immunohistochemistry	Conforms
----------------------	----------

working dilution 1:50

Specification: PRD.1.ZQ5.10000054505

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.