

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Linoleic acid - ≥99%

Product Number: L1376
CAS Number: 60-33-3
MDL: MFCD00064241
Formula: C18H32O2
Formula Weight: 280.45 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Purity (GC)	≥ 99 %
Proton NMR spectrum	Conforms to Structure
Recommended Retest Period	-----
4 years	

Specification: PRD.0.ZQ5.10000034933

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.