

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lyticase from *Arthrobacter luteus* – lyophilized powder, $\geq 2,000$ units/mg protein, Protein ≥ 20 % by biuret

Product Number:

L2524

CAS Number:

37340-57-1

MDL:

MFCD00131560

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

% Protein (Biuret)

 ≥ 20

units/mg protein

 ≥ 2000

One unit will produce a change in absorbance at 800nm of 0.001 per minute at pH 7.5 at 25 deg C, using a suspension of brewer's yeast as substrate in a 3 ml reaction mixture.

Recommended Retest Period

1 year

Specification: PRD.0.ZQ5.10000030129

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.