

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(ethylene glycol) linear dendron - NHBoc terminated, PEG Mn 2000, generation 4

Product Number:

901405

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Viscous Liquid

1H NMR Spectrum

Conforms to Structure

Molecular Number by GPC

3500 - 8500

Recommended Retest Period

2 Years

Specification: PRD.0.ZQ5.10000101988

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.