

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Pro-Phe-Arg-7-amido-4-methylcoumarin - ≥95%

**Product Number:****P9273**

CAS Number:

65147-21-9

MDL:

MFCD00038452

Formula:


C30H37N7O5

Formula Weight:

575.66 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

White to Off White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear

50 mg/mL, MeOH

Peptide Content

≥ 65 %

Purity (HPLC)

≥ 95 %

Specification: PRD.1.ZQ5.10000035595

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.