

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-phospho-Synapsin I (pSer<sup>9</sup>) antibody produced in rabbit – ~1 mg/mL, affinity isolated antibody, buffered aqueous solutionProduct Number: **S8067**

MDL: MFCD03792020

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in 10 mM Hepes, pH 7.5, containing 150 mM NaCl, 100 ug/mL BSA, 50% Glycerol	
Working Dilution	Conforms
By Direct Immunoblotting; 1:1000	

Specification: PRD.0.ZQ5.10000090176

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.