


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
6-Hydroxypyridine-2-carboxylic acid - 95%

Product Number: 384305
CAS Number: 19621-92-2
MDL: MFCD00192220
Formula: C₆H₅NO₃
Formula Weight: 139.11 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Titration with NaOH

Very Faint Yellow to Yellow and Very Faint Brown to Brown
Powder or Crystals
Conforms to Structure
94.5 - 105.5 %

Specification: PRD.2.ZQ5.10000029450

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.