

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Linolenic acid - $\geq 99\%$

Product Number: L2376
CAS Number: 463-40-1
MDL: MFCD00065720
Formula: C18H30O2
Formula Weight: 278.43 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Purity (GC)	$\geq 99\%$
Proton NMR spectrum	Conforms to Structure
Recommended Retest Period	-----
4 years	

Specification: PRD.2.ZQ5.10000040125

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.