

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2-Bis(2-Aminophenoxy)ethane-N,N,N',N'-tetraacetic acid - 98%

Product Number: A4926
CAS Number: 85233-19-8
MDL: MFCD00036255
Formula: C₂₂H₂₄N₂O₁₀
Formula Weight: 476.43 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow to Faint Red
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/ml, 0.3 N NaHCO ₃	
Carbon	54.4 - 56.6 %
Nitrogen	5.7 - 6.0 %
¹ H NMR Spectrum	Conforms to Structure
¹³ C NMR Spectrum	Conforms to Structure
Recommended Retest Period	-----
3 years	

Specification: PRD.1.ZQ5.10000007508

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.