


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
tert-Butyl diazoacetate

Product Number: 480754
CAS Number: 35059-50-8
MDL: MFCD01318409
Formula: C₆H₁₀N₂O₂
Formula Weight: 142.16 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Clarity)	Clear to Hazy
Appearance (Color)	Colorless to Yellow and Colorless to Orange
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 80.0 %
Dichloromethane	≤ 15.0 %

Specification: PRD.0.ZQ5.10000055388

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.