


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Iodoacetic anhydride

Product Number: 284262
CAS Number: 54907-61-8
MDL: MFCD00001080
Formula: C₄H₄I₂O₃
Formula Weight: 353.88 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Dark Yellow
Appearance (Form)	Conforms
Powder, Flakes, Crystals and/or Chunks	
Solubility (Turbidity)	Clear to Slightly Hazy
C = 10%, CHCl ₃	
Solubility (Color)	Conforms
Yellow, Pink or Tan	
Proton NMR Spectrum	Conforms to Structure

Specification: PRD.1.ZQ5.10000037390

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.