


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Antimony(III) sulfide - powder

Product Number: 244562
CAS Number: 1345-04-6
MDL: MFCD00011217
Formula: S3Sb2
Formula Weight: 339.72 g/mol



TEST	Specification
Appearance (Color)	Grey to Very Dark Grey
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Titration with Iodine % Sb	≥ 68.5 %

Specification: PRD.1.ZQ5.10000003798

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.