


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Arg-Gly-Asp-Ser-Pro-Ala-Ser-Ser-Lys-Pro - $\geq 97\%$ (HPLC)

Product Number: A6677
CAS Number: 91575-25-6
MDL: MFCD00076453
Formula: C₄₀H₆₈N₁₄O₁₆
Formula Weight: 1,001.05 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/ml, H ₂ O	
Peptide Content	$\geq 65\%$
Purity (HPLC)	$\geq 97\%$

Specification: PRD.0.ZQ5.10000054973

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.