

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Biotin disulfide N-hydroxysuccinimide ester - ≥90% (TLC), powder

Product Number:**B4531**

CAS Number:

142439-92-7

MDL:

MFCD00466919

Formula:


C₂₂H₃₃N₅O₇S₃

Formula Weight:

575.72 g/mol

Storage Temperature:

-20 °C



TEST	Specification
------	---------------

Purity (HPLC)

≥ 98 %

Mass Spectrum

Conforms

Consistent with Structure

Specification: PRD.1.ZQ5.10000014781

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.