


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lipoxin A4 - ethanol solution

Product Number: **L0521**
CAS Number: 89663-86-5
MDL: MFCD00065845
Formula: C₂₀H₃₂O₅
Formula Weight: 352.47 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Concentration by UV-VIS in EtOH	49 - 350 µg/ml
Purity (HPLC)	≥ 92.5 %

Specification: PRD.1.ZQ5.10000060199

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.