

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Vinyl acetate - contains 3-20 ppm hydroquinone as inhibitor, ≥99%

Product Number:

V1503

CAS Number:

108-05-4

MDL:

MFCD00008713

Formula:


C₄H₆O₂

Formula Weight:

86.09 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 99.0 %

HQ as Inhibitor

3 - 20 ppm

Specification: PRD.0.ZQ5.10000017301

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.