

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Glucosidase from almonds - lyophilized powder, 10-30 units/mg solid

Product Number:

G4511

CAS Number:

9001-22-3

TEST	Specification
------	---------------

Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White to Light Yellow/Tan
Solubility (Solvent)	Water
Solubility (Conc)	1.90 - 2.10 mg/ml
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless to Faint Yellow
Protein Content (Method)	Biuret
Prot. Content (% Protein)	≥ 90.00 %
Enzym. Assay (Activity U/mgS)	10.00 - 30.00 U/mgS
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000001384

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.