

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Aminopeptidase I from Streptomyces griseus - lyophilized powder,  $\geq 200$  units/mg protein

Product Number:

**A9934**

CAS Number:

9031-94-1

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White to Faint Yellow
Protein Content (Method)	Lowry
Prot. Content (% Protein)	40.00 - 60.00 %
Enzym. Assay (Spec. Activity)	$\geq 200.00$ U/mgP
Vial Content (Size)	100.00 - 100.00
Vial Content (Unit Of Measure)	$\mu\text{g}$ Protein
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000002015

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.