

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Aluminum acetylacetonate ReagentPlus®, 99%

Product Number:

208248

CAS Number:

13963-57-0

MDL:


MFCD00000013

Formula:

C₁₅H₂₁AlO₆

Formula Weight:

324.31 g/mol

**TEST****Specification**

Appearance (Color)

White to Yellow

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Complexometric EDTA

8.1- 8.5 %

% Al

Remarks:

Specification Date : 12/16/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.