## SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com


## **Product Specification**

AU-s O- OH

Product Name: Gold nanoparticles – 5 nm diameter, carboxylic acid functionalized, PEG 3000 coated, OD 50, dispersion in H2 O

**Product Number:** 

765430



TEST	Specification	
Appearance (Color)	Red to Very Dark Red	
Appearance (Form)	Liquid	
Precipitate Present	None Present	
Agglomeration/Precipitate		
Identity	Confirmed	
Gold Nanoparticles, Carboxylic Acid Functionalized, PEG 3000 Coated		
Dispersion in Water	Confirmed	
Concentration	Conforms	
(OD) = 50		
Core Size	3 - 7 nm	
(Diameter)		
Recommended Retest Period		
12 Months		

Specification: PRD.0.ZQ5.10000060112

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.