

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Renin Substrate 1

Product Number: R8276

MDL: MFCD00284486

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	Orange to Dark Red
Appearance (Form)	Powder
Solubility (Color)	Orange to Red
Solubility (Turbidity)	Clear
1 mg/mL, H ₂ O	
Peptide Content	≥ 60 %
BY AAA	
Purity (HPLC)	≥ 85 %

Specification: PRD.0.ZQ5.10000059882

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.