


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2-Oxopiperazine - 97%

**Product Number:** 641065  
**CAS Number:** 5625-67-2  
**MDL:** MFCD01318687  
**Formula:** C<sub>4</sub>H<sub>8</sub>N<sub>2</sub>O  
**Formula Weight:** 100.12 g/mol



### TEST

### Specification

Appearance (Color)  
Appearance (Form)  
Powder, Crystals or Crystalline Powder and/or Chunks  
Infrared spectrum  
Carbon  
Nitrogen  
Purity (GC)

Faint Yellow to Dark Yellow and Faint Orange to Dark Orange  
Conforms to Requirements

Conforms to Structure  
46.3 - 49.7 %  
27.0 - 29.0 %  
≥ 96.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.