

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Betulinic acid - ≥98% (HPLC)


Product Number: **B8936**

CAS Number: 472-15-1

Formula: C₃₀H₄₈O₃

Formula Weight: 456.70 g/mol

Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/mL, EtOH	
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 2 %
Purity (HPLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000039978

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.