

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Decaethylene glycol mono-dodecyl ether - nonionic surfactant

Product Number:

P9769



CAS Number:

9002-92-0

MDL:

MFCD00080891

Formula:

C32H66O11

Formula Weight:

626.86 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Semi-Solid
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, EtOH	
Acid Value	≤ 0.15
Hydroxyl Value	85 - 95

Specification: PRD.0.ZQ5.10000036828

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.