

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Bis(cyclopentadienyl)cobalt(III) hexafluorophosphate 98%

Product Number:

279811

CAS Number:

12427-42-8

MDL:


MFCD00011874

Formula:

C₁₀H₁₀CoF₆P

Formula Weight:

334.08 g/mol



TEST

Specification

Appearance (Color)
Yellow to Orange to Brown
Appearance (Form)
Powder or Crystals
Infrared spectrum

Carbon

Conforms to Requirements
Conforms to Requirements
Conforms to Structure

35.0- 36.8 %

Remarks:

Specification Date : 11/13/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.