


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
1,3-Diamino-2-propanol - 95%

**Product Number:** D18609  
**CAS Number:** 616-29-5  
**MDL:** MFCD00008142  
**Formula:** C<sub>3</sub>H<sub>10</sub>N<sub>2</sub>O  
**Formula Weight:** 90.12 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Colorless/White to Yellow to Beige	
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Purity (Titration by HClO <sub>4</sub> )	94.0 - 106.0 %
Purity (GC)	≥ 94.0 %
Solubility (Turbidity)	Clear to Slightly Hazy
c = 50 mg/ml (5%); Water	
Solubility (Color)	Colorless to Light Yellow

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.