


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Nafion® perfluorinated resin, powder

Product Number: **495786**  
CAS Number: 31175-20-9  
MDL: MFCD00132832



TEST	Specification
Appearance (Color)	White to Tan
Appearance (Form)	Crystalline Powder
Equivalent Weight g/Eq, with NaOH	990 - 1300
Registered Trademark	Confirmed
Nafion® is a registered trademark of The Chemours Company FC LLC	

Specification: PRD.0.ZQ5.10000061379

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.