

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lysostaphin from *Staphylococcus staphylolyticus* - lyophilized powder, Protein 50-70 % by biuret,
≥500 units/mg protein

Product Number:

L7386

CAS Number:

9011-93-2

TEST**Specification**

Appearance (Form)	Lyophilized Powder
Appearance (Colour)	White to Light Yellow
Solubility (Solvent)	Water
Solubility (Conc)	2.00 - 2.10 mg/ml
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless to Light Yellow
Protein Content (Method)	Biuret
Prot. Content (% Protein)	50.00 - 70.00 %
Enzym. Assay (Spec. Activity)	≥ 500.00 U/mgP
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000001223

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.