

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Tris-EDTA buffer solution - 100X concentrate for Northern and Southern blotting

Product Number:

T9285



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
pH	7.8 - 8.2
1x, H ₂ O	
DNAse, Exonuclease Detection	None Detected
RNAse Detection	None Detected
Conductivity	≥ 30000
Measured in umhos	
Standard Plate Count	≤ 10 CFU/ml
Recommended Retest Period	-----
3 years	

Specification: PRD.2.ZQ5.10000006975

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.