

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

13-cis-Retinoic acid - ≥98% (HPLC)

 Product Number:
 R3255

 CAS Number:
 4759-48-2

 MDL:
 MFCD00079542

 Formula:
 C20H28O2

 Formula Weight:
 300.44 g/mol

Storage Temperature: -20 ℃

GC CH3 CH5 CH5

TEST

Appearance (Color)
Appearance (Form)
Infrared spectrum
Carbon

Carbon
Purity (HPLC)
Solubility (Color)
Solubility (Turbidity)
50 mg/mL, CHCl3
Recommended Retest Period

3 Years

Specification: PRD.1.ZQ5.10000043200

Specification

Yellow and Orange and Yellow-Orange

Pow der

Conforms to Structure

78.4 - 81.6 %

≥ 98 %

Yellow to Yellow-Green

Clear

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.