


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Methoxy-TEMPO - for ESR-spectroscopy

Product Number: 514586
CAS Number: 95407-69-5
MDL: MFCD00270334
Formula: C₁₀H₂₀NO₂
Formula Weight: 186.27 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Orange to Red
Appearance (Form)	Conforms to Requirements
Solid, Powder or Crystals	
Infrared spectrum	Conforms to Structure
Carbon	62.2 - 66.7 %
Nitrogen	7.2 - 7.8 %
Purity (HPLC)	≥ 96.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.