

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:
3-Hydroxypropionic acid, 30% solution in water

Product Number: 792659

Formula: C₃H₆O₃

Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Proton NMR spectra	Conforms to Structure
Water (by Karl Fischer)	65.0 - 75.0 %
Titration with NaOH (anhydrous basis)	84.5 - 115.5 %

Specification: PRD.0.ZQ5.10000059805

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.