

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Zinc perchlorate hexahydrate

Product Number: 401439
CAS Number: 10025-64-6
MDL: MFCD00150370
Formula: Cl₂O₈Zn · 6H₂O
Formula Weight: 372.38 g/mol

Zn(ClO₄)₂ · 6H₂O

TEST	Specification
Appearance (Color)	White
Appearance (Form) Crystals or Crystalline Powder and/or Chunks	Conforms to Requirements
X-Ray Diffraction	Conforms to Structure
Complexometric EDTA % Zinc	

Specification: PRD.1.ZQ5.10000101684

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.